


## X-CLARITY™ Polymerization Initiator C13104



### Storage

4°C

### Product Description

Molecular formula	C <sub>12</sub> H <sub>22</sub> N <sub>4</sub> · 2HCl
Molecular weight	323.33 g/mol
Appearance	Pale yellow powder
10-hr half-life	44°C (in water)

The X-CLARITY™ Polymerization Initiator is a thermal free radical initiator also known as VA-044. VA-044 is relatively stable at room temperature but undergoes rapid hemolytic decomposition when heated in aqueous solution, releasing cationic free radicals that initiate the polymerization of hydrogel monomers.

### Directions for Use

#### STOCK PREPARATION

1. Dissolve 2.5 g X-CLARITY™ Polymerization Initiator in 10 mL 1X PBS to make a 25% (w/v) stock solution.
2. Aliquot and store at -80°C for up to 6 months.
3. Thaw at 4°C or on ice before use.

#### HYDROGEL INFUSION

1. Prepare fresh hydrogel-initiator solution as needed.
  - 4% acrylamide (A4P0): Add one part 25% (w/v) X-CLARITY™ Polymerization Initiator to 100 parts X-CLARITY™ Hydrogel Solution. Mix thoroughly.
  - 1% acrylamide (A1P0): Make the 4% acrylamide solution above and dilute 1:4 in 1X PBS. Mix thoroughly.
2. Use enough mixture to fully submerge samples.
3. Incubate at 4°C for 24 hours.

#### TISSUE-HYDROGEL HYBRIDIZATION

1. Initiate polymerization with the X-CLARITY™ Polymerization System. Run the system at 37°C for 3 hours at -90 kPa.
2. Shake samples gently on a shaker for 1 minute. If in a conical tube, invert the sample gently for 1 minute.

! *Bis*-acrylamide creates crosslinks between polyacrylamide chains, which hardens the hydrogel network and forms a rigid gel around tissue samples that must be removed prior to clearing. X-CLARITY™ Hydrogel Solution does not contain *bis*-acrylamide, which prevents a gel from forming around the sample. This is ideal for small or delicate tissues. A successfully polymerized *bis*-acrylamide-free hydrogel solution is sticky. The final shaking step is important to ensure a homogenous distribution of the solution throughout the sample.

### Disclaimer

This product is for research use only. Please consult the material safety data sheet for information regarding hazards and safe handling practices.

### References

1. Lee, E et al. ACT-PRESTO: Rapid and consistent tissue clearing and labeling method for 3 dimensional (3D) imaging. *Sci Rep* 6, 18631 (2016).
2. Yang, B. et al. Single-cell phenotyping within transparent intact tissue through whole-body clearing. *Cell* 158, 945–958 (2014).
3. Chung, K. et al. Structural and molecular interrogation of intact biological systems. *Nature* 497, 332–337 (2013).

Additional information is available on our website at [www.logosbio.com](http://www.logosbio.com). © Logos Biosystems 2019. All rights reserved.



#### HEADQUARTERS

FL 2 & 3  
28 Simindaero 327beon-gil, Dongan-gu  
Anyang-si, Gyeonggi-do 14055  
South Korea

Tel: +82 (31) 478-4185

#### USA

7700 Little River Turnpike STE 207  
Annandale, VA 22003  
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

#### EUROPE

11B avenue de l'Harmonie  
59650 Villeneuve d'Ascq  
France

Tel: +33 (0)3 74 09 44 35