

# Yeast Cell Counting using the LUNA-FX7™ Automated Cell Counter

**Key features:** Accurate yeast cell counts using the LUNA-FX7™ Automated Cell Counter

## INTRODUCTION

Yeasts are used in a vast array of processes including basic research, brewing & distilling, and food production that require accurate cell counts and viabilities prior to and throughout the process.<sup>1,2,3</sup>

The small physical size, small genome size, morphology, and development have proven to be quite challenging for automated cell counters to overcome.<sup>4,5</sup> As a result, the tedious and error-prone method of manually counting yeast using a hemocytometer and vital stain remains the 'go to' method. Here we show how the LUNA-FX7™ can accurately and reliably count yeast cells.

## MATERIALS

1	LUNA-FX7™ Automated Cell Counter	5	LUNA™ Cell Counting Slide [ L12001 ] / PhotonSlide™ [ L12005 ]
2	EC1118	6	Cell Dilution Buffer II [ F53002 ]
3	0.02 % Methylene Blue [ L13003 ]	7	Hemocytometer
4	Acridine Orange [ F23002 ] / Propidium Iodide [ F23003 ]		

## METHODS

- 1) Yeast strain EC1118 was inoculated in YPD broth and incubated overnight at 30 °C, 200 rpm
- 2) Yeast was diluted in 1XPBS at ratios of 1:3 to 1:5.
- 3) Diluted cultures were diluted into 2 equal samples. One sample was then heat killed (65 °C for 20 minutes) to generate 'Live and Dead' cells for viability testing.
- 4) Cells were further diluted using Logos Biosystems Yeast Cell Dilution Buffer II : 1:5 and 1:10 (cells: buffer)

### Viability Testing

- 5) Live and dead cells were mixed as per the Table1

**Table 1.** Mixture ratio of live cells and dead cells for viability test

	0%	25 %	50 %	75 %	100 %
Live Cells	0 µl	25 µl	50 µl	75 µl	100 µl
Dead Cells	100 µl	75 µl	50 µl	25 µl	0 µl


## RESULTS

### **Linearity of viability of the LUNA-FX7™ Automated Cell Counter**


To determine the linearity of viability, samples of various viability were prepared using live cells and dead cells and then counted using the LUNA-FX7™ Automated Cell Counter. With the LUNA-FX7™, viability can be checked in fluorescence mode and [CELL COUNTING & VIABILITY] of brightfield mode.

To examine viability in brightfield, yeast samples were mixed with 0.02 % methylene blue stain in the same volume, and 10ul was loaded on a LUNA™ Cell Counting Slide. Yeast samples were prepared as dead cells (0 %), 25 %, 50 %, 75 %, live cells (100 %). Referring to Figure 1A, five samples were stained with 0.02 % methylene blue stain, and images taken in the [CELL COUNTING & VIABILITY] of brightfield can be observed. Counting was repeated 3 times to calculate average viability, and the linearity of viability was confirmed accordingly. As shown in Figure 1B, viability measured by [CELL COUNTING & VIABILITY] of brightfield mode of LUNA-FX7™ was  $R^2=0.9964$ , which showed high agreement with the theoretical survival rate.

**A**




**B**




**Figure 1.** Cell viability determined with [CELL COUNTING & VIABILITY] of brightfield mode of LUNA-FX7™ Automated Cell Counter. (A) It shows images of yeast cells stained with methylene blue and tagged. (B) The measured viability and the theoretically expected viability were confirmed with linear graph. ( $R^2=0.9964$ )

In addition, to investigate viability in fluorescence, five yeast samples were stained with AO/PI, and 10ul was loaded onto a PhotonSlide™. [CELL LINES & PRIMARY CELLS] had been the only option in the fluorescence mode of LUNA-FX7™. Through the last software update, a new fluorescence mode was added [CELL LINES & PRIMARY CELLS, ADVANCED] to provide better results for cells with a lot of aggregation such as yeast. After AO/PI staining, images were acquired and observed as shown in Figure 2A. As seen in Figure 2B, five different viability samples were counted 3 times, and even in [ADVANCED] mode, the measured viability of yeast cells was  $R^2=0.9979$ , which had a significant correlation with the theoretical viability.

**A**



**B**




**Figure 2.** Cell viability determined with [CELL LINES & PRIMARY CELLS, ADVANCED] of fluorescence mode of LUNA-FX7™ Automated Cell Counter. (A) It shows images of yeast cells stained with acridine orange/propidium iodide and tagged. (B) The measured viability and the theoretically expected viability were confirmed with linear graph. ( $R^2=0.9979$ )

### **Linearity of concentration analysis of the LUNA-FX7™ Automated Cell Counter**


To determine the linearity at various concentration in counting yeast using LUNA-FX7™ Automated Cell Counter, high concentration yeast was prepared, and the concentration was measured by serial dilution. For observation in

brightfield mode, it was stained with 0.02 % methylene blue, and counted in [TOTAL CELL COUNTING] and [CELL COUNTING & VIABILITY] in brightfield mode, respectively. As can be seen in Figure 3, the  $R^2$  value were measured to be 0.998 and 0.9991, respectively. This is similarly observed in the fluorescence mode, and when it was stained with AO/PI and captured in [ADVANCED] mode, as shown in figure 4, the  $R^2$  value was 0.9753. This means that yeast counting using LUNA-FX7™ shows high linearity from high to low concentrations, and it means that more accurate measurement is possible with LUNA-FX7™.

Through these experiments and data from the previous application note [An analysis of optimal cell counting range of the LUNA-FX7™], we can think about the limit of quantification about yeast cells. The average size of yeast cells is 3~4  $\mu\text{m}$ , and the size of a typical mammalian cell is 10  $\mu\text{m}$ . Therefore, compared to the range of  $5.00 \times 10^4$  to  $1.50 \times 10^7$  cells/mL for mammalian cells in the 2-ch slide, it can be adjusted to  $5.00 \times 10^5$  to  $3.00 \times 10^7$  cells/mL for yeast cells. Because the cell size differs by more than 2 times, it is impossible to count at a low concentration of  $5.00 \times 10^4$  cells/mL as in mammalian cells, but due to the size difference, it is possible to count cells with a high concentration that is twice as high as  $1.50 \times 10^7$  cells/mL.




**Figure 3.** Concentration of yeast sample measured in brightfield mode of LUNA-FX7™ Automated Cell Counter. (A) [TOTAL CELL COUNTING] (B) [CELL COUNTING & VIABILITY]



**Figure 4.** Concentration of yeast sample measured in [CELL LINES & PRIMARY CELLS, ADVANCED] of fluorescence mode of LUNA-FX7™ Automated Cell Counter.

### **LUNA-FX7™ Automated Cell Counter vs. Hemocytometer**

An experiment was conducted to compare the LUNA-FX7™ and hemocytometer, which have been traditionally used. As in the previous experiment, the sample was stained with 0.02 % methylene blue and AO/PI, and 10 µl each was loaded on a hemocytometer, LUNA™ Cell Counting Slide and PhotonSlide™. The images were taken in two modes of brightfield and fluorescence mode of LUNA-FX7™, and the concentrations were compared with those measured on a hemocytometer. As shown in Figure 5, no significant difference was observed between the values measured by the hemocytometer and the values measured by the LUNA-FX7™. When comparing the values in the hemocytometer and each mode, the P values are all over 0.05, and there is no statistically significant difference.



**Figure 5.** Comparison of yeast counting results of hemocytometer and LUNA-FX7™ Automated Cell Counter.

### **CONCLUSION**

Here, we've shown that the LUNA-FX7™ can accurately and reliably count yeast cells. Further when compared to the traditional manual method, counting yeast with the LUNA-FX7™ is 20X faster, more reliable, and less prone to error.

## REFERENCES

- <sup>1</sup> Legras JL, Merdinoglu D, Cornuet JM, Karst F (2007). "Bread, beer and wine: *Saccharomyces cerevisiae* diversity reflects human history". *Molecular Ecology*. 16 (10): 2091–2102. doi:10.1111/j.1365-294X.2007.03266.x. PMID 17498234. S2CID 13157807.
- <sup>2</sup> "Fuel Ethanol Production: GSP Systems Biology Research". Genomic Science Program. U.S. Department of Energy Office of Science. Archived from the original on 3 June 2009. Retrieved 28 November 2009.
- <sup>3</sup> Barnett JA (2003). "Beginnings of microbiology and biochemistry: the contribution of yeast research" (PDF). *Microbiology*. 149 (3): 557–567. doi:10.1099/mic.0.26089-0. PMID 12634325. S2CID 15986927. Archived from the original (PDF) on 3 March 2019.
- <sup>4</sup> Walker K, Skelton H, Smith K (2002). "Cutaneous lesions showing giant yeast forms of *Blastomyces dermatitidis*". *Journal of Cutaneous Pathology*. 29 (10): 616–618. doi:10.1034/j.1600-0560.2002.291009.x
- <sup>5</sup> Britannica, The Editors of Encyclopedia. "yeast". *Encyclopedia Britannica*, 13 May. 2020, <https://www.britannica.com/science/yeast-fungus>. Accessed 22 December 2021.