

FTP FREEWARE – FileZilla

PC Requirements

OS: Windows 7, Windows 10
Display: 1280 x 1024 pixels
Internet: Wi-Fi Connection



HEADQUARTERS

FL 2 & 3 28
Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055

SOUTH KOREA

TEL: +82 (31) 478-4185

7700 Little Turnpike Street
Suite 207

Annandale, VA 22003

UNITED STATES

TEL: +1 (703) 622-4660

11B avenue de l'Harmonie
59650 Villeneuve d'Ascq

FRANCE

TEL: +33 (0)3 7 4 09 44 35

www.logosbio.com

REV3_20AUG21

Product Description

FTP Freeware, *FileZilla*, is software used to access data saved to the QUANTOM Tx™ Microbial Cell Counter from a remote PC via Wi-Fi.

Directions for Use

CONNECT THE QUANTOM Tx™ TO Wi-Fi

1. Connect the supplied Wi-Fi dongle to a USB port.
2. Open Wi-Fi Options. Log in to a wireless network.
3. Make a note of the QUANTOM Tx™ IP address.

INSTALL FTP FREEWARE, *FileZilla* ON A PC

1. Locate the serial number of the QUANTOM Tx™

The serial number can be found on the back of the unit. It will begin with the letters 'QTX' followed by seven numbers. (e.g.: QTX-00-00000).

2. QUANTOM Tx™ units with serial numbers ending in -00039 or higher will proceed with the installation of *FileZilla*.

! If the serial number of the QUANTOM Tx™ ends in -00001 through -00038, please refer to the QUANTOM Tx™ User Manual or the technote on 'Using QUANTOM™ Utility to transfer data via Wi-Fi for QUANTOM Tx™'!


3. Download the FTP Freeware, *FileZilla* from filezilla-project.org.
4. Unzip the downloaded file and install on a compatible PC.

CONNECT *FileZilla* TO THE QUANTOM Tx™

1. Download and install **FileZilla Client** from filezilla-project.org.
2. Open **FileZilla Client**. To connect the QUANTOM Tx™ enter the following:
 - Host: Enter the QUANTOM Tx™ IP address
 - Username: logos
 - Password: logos

! The QUANTOM Tx™ and the PC must be on the same wireless network!

3. Enter the QUANTOM Tx™ IP address. Click **QuickConnect**.



4. The **REMOTE SITE** section will display the folders/files saved to the QUANTOM Tx™.
5. Select the desired folders/files and click **DOWNLOAD** to import to the PC.

For more detailed instruction refer to the product manual.
© Logos Biosystems 2021. All Rights Reserved.