













Vessel holders

A collection of vessel holders have been to accommodate various flasks, dishes, plates, and slides.

For use with the Vessel Holder Frame [Cat # I10230 – I10241]

<p>I10230 Universal vessel holder</p>  <p>Ø30</p>	<p>I10231 Microplate vessel holder</p>  <p>Ø128.3 Ø86</p>	<p>I10232 Single slide vessel holder</p>  <p>75 25</p>
<p>I10233 35 mm Cell Culture Dish Holder</p>  <p>Ø35</p>	<p>I10234 60 mm Cell Culture Dish Holder</p>  <p>Ø60</p>	<p>I10235 100 mm Cell Culture Dish Holder</p>  <p>Ø100</p>
<p>I10237 Nunc T-25 Flask Holder</p>  <p>Ø77.5 Ø47.5</p>	<p>I10238 Nunc T-75 Flask Holder</p>  <p>Ø128.5 Ø79.5</p>	<p>I10239 Falcon T-25 Flask Holder</p>  <p>Ø93 Ø41</p>
<p>I10240 Falcon T-75 Flask Holder</p>  <p>Ø143 Ø86</p>	<p>I10241 Nunc CellCulture Plate Holder</p>  <p>Ø128.3 Ø86</p>	<p>I10242 PerkinElmer Microplate Holder</p>  <p>Ø128.3 Ø86</p>



www.logosbio.com

HEADQUARTERS

FL 2 & 3
28 Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055
South Korea

Tel : +82 (31) 478-4185
Fax : +82 (31) 360-4277
E-mail : info@logosbio.com

USA

7700 Little River Turnpike STE 207
Annandale, VA 22003
USA

Tel : +1 (703) 622-4660, +1 (703) 942-8867
Fax : +1 (571) 266-3925
E-mail : info-usa@logosbio.com

EUROPE

118 avenue de l'Harmonie
59650 Villeneuve d'Ascq
France

Tel : +33 (0)3 74 09 44 35
Fax : +33 (0)3 59 35 01 98
E-mail : info-france@logosbio.com