

LUNA™ Standard Beads

B13101

Storage

4 °C

Product Description

Concentration	~1.0 x 10 ⁶ beads/mL (at dilution factor 2)
Average size	9-11 µm
Appearance	Colorless beads in dark blue liquid
Bead content	Polystyrene
Dye content	Trypan blue

LUNA™ Standard Beads are used to calibrate the LUNA™ family of cell counters for brightfield cell counting.

Directions for Use

Mix LUNA™ Standard Beads thoroughly before use to ensure a homogenous suspension. Work quickly when using LUNA™ Standard Beads to ensure the beads do not settle.

Disclaimer

This product is for research use only.
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website at www.logosbio.com.
© Logos Biosystems 2018. All rights reserved.



HEADQUARTERS

FL 2 & 3
28 Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055
South Korea

Tel: +82 (31) 478-4185

USA

7700 Little River Turnpike STE 207
Annandale, VA 22003
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

EUROPE

118 avenue de l'Harmonie
59650 Villeneuve d'Ascq
France

Tel: +33 (0)3 74 09 44 35

www.logosbio.com
VL1801-01