

CountWire™

Version Release Note: 1.2.0

CountWire



Dec. 9, 2021

Sales & Marketing Department

< Version 1.2.0 >

Release Date	Model	Version	Writer
Dec. 9, 2021	CountWire™	1.2.0	Deokhui Chang

Hardware	Software	New features	Improvement
-	o	o	o

Version description

1. Software bugs

- A bug has been fixed where the system failed to automatically log out the LUNA-FX7 user when password was changed or reset in CountWire. In the revised version, the users are automatically logged out.
- A bug has been fixed in which Printed names longer than 16 characters were truncated.

2. Modifications

- The software downgrade of Data Storage prohibited due to possible inconsistency and complexity.

3. Improvements

- (Administrator menu only) File transfer option has been added in **System**. For "**Transfer PDF reports only to save storage space**", only PDF reports and the first field images are transferred to Data Storage to save storage space. For "**Transfer PDF reports and all image files**", both PDF reports and image files are transferred to Data Storage.
- (Administrator menu only) New back up current storage information & register new storage information menu has been added in the System menu. By pressing START of **Backup current storage information**, it stores the current data storage information and by pressing START of **Registrar new storage information**, it restores the previous data storage information to the new hard drive.
- (Administrator menu only) A remote access feature has been added for future technical support.
- **Staining option** has been removed from **Instrument setting > Cell counting setting**.
- The IQOQ protocol has been added for all counting modes.
- The YEAST protocol has been added for Brightfield and Fluorescence Advanced modes.

- <Caps Lock> indication has been added in login window and in **Creat account..**

The image displays two side-by-side screenshots of the CountWire software interface. The left window is the login screen, featuring the CountWire logo on the left and a form on the right. The form includes fields for 'Storage IP' (with sub-fields 192, 168, 0, 176), 'Storage port' (22), 'User ID' (administrator), and 'Password' (masked with dots). A 'Caps Lock is On' indicator is visible below the password field. There are 'LOG IN' and 'SIGN UP' buttons at the bottom. The right window is the 'Create account' screen, with fields for 'User ID', 'Print name', 'Password *', and 'Confirm *'. The 'Password *' field has a 'Caps Lock is On' indicator and a note '* At least 8 characters'. There are 'OK' and 'CANCEL' buttons at the bottom.

- End of Page -