

LUNA-FX7™ Automated Cell Counter

Version Release Note: 1.9.7



Sep. 25, 2025

Sales & Marketing Department

< Version 1.9.7 >

Release Date	Model	Version	Writer
Sep 25, 2025	LUNA-FX7™	1.9.7	Yubeen Lee

Hardware	Software	New features	Improvement
-	0	-	0

Version description

<Improvement>

1. WiFi – Static IP

- In versions 1.9.2 to 1.9.6, the Wi-Fi Static IP option was disabled and unavailable.
- This option has now been re-enabled and is fully functional again.

<Bug fix>

1. Security mode – Physical Keyboard Shift Input

- Previously, when using a physical keyboard during login, pressing the Shift key together with another key could cause the Shift key itself to be entered as a character.
- This issue has been fixed so that only the intended special character is correctly entered.

<Note for CountWire customers!>

- To ensure compatibility, please update **CountWire Client** and **DSS** to the versions specified in the table when upgrading **LUNA-FX7** to this version.

Component	Required Version
CountWire Client	v1.3.6
CountWire DSS	v1.3.5

Download the latest software via website <https://logosbio.com/support/software>

- End of Page -