

LUNA-FX7™ Automated Cell Counter

Version Release Note: 1.9.1



Jun. 30, 2023

Sales & Marketing Department


< Version 1.9.1 >

Release Date	Model	Version	Writer
Jun. 30, 2023	LUNA-FX7™	1.9.1	Yubeen.Lee
Hardware	Software	New features	Improvement
-	0	-	0

Version description

<Improvement>

1. Autofocus algorithm optimization
2. Add Autofocus Offset setting for user convenience



*Download the latest software via website <https://logosbio.com/support/software>

- End of Page -