

LUNA-III™ Automated Cell Counter

Version Release Note: 1.0.1



Sales & Marketing Department

<LUNA-III™ Version Release History>

Version	Release date	Version Summary	Page
1. 0. 1	Dec. 16, 2024	· The bugs in live view and software update have been fixed.	1

< Version 1.0.1 >

Release Date	Model	Version	Writer
16-Dec-2024	LUNA-III™	1.0.1	Yubeen Lee

Hardware	Software	New features	Improvement
-	0	-	-

Version description
<ol style="list-style-type: none">1. Fixed S/W bug<ul style="list-style-type: none">- Fixed a bug where image update was slow during live view.- Fixed a bug where the software update date and time were incorrect.- Fixed a bug where the scalebar would sometimes not displayed on tagged image.- Fixed a bug where the focus position would gradually increase during focus calibration.

- End of Page -