

# QUANTOM Tx™ Data Transfer Utility Installation manual v1.0









## **IMPORTANT**

- We support two Data Transfer Utility to import QUANTOM's data to PC.
1. **QUANTOM Tx Utility** is only available for Windows.  
(To use it, install the program located in the QUANTOM Tx Utility folder.)
  2. **FTP freeware (FileZilla Client)** has Windows, Mac, and Linux versions, but it is recommended to use it in a Windows environment. (You can download the installation file from the website <https://filezilla-project.org/>)

# QUANTOM Tx™ Data Transfer Utility Installation manual


## ❖ QUANTOM Tx™ Utility Installation manual

- Check your system information and install the version that matches the bit of your operating system. Use **QUANTOM-Utility-win64-v1.2.0.exe** if your PC System type is a **64-bit Operating System** and **QUANTOM-Utility-win32-v1.2.0.exe** if your PC System type is a **32-bit Operating System**.



# QUANTOM Tx™ Data Transfer Utility Installation manual

## ❖ FTP freeware(FileZilla Client) Installation manual



Please select your edition of FileZilla Client

	FileZilla	FileZilla with manual	FileZilla Pro	FileZilla Pro + CLI
Standard FTP	Yes	Yes	Yes	Yes
FTP over TLS	Yes	Yes	Yes	Yes
SFTP	Yes	Yes	Yes	Yes
Comprehensive PDF manual	-	Yes	Yes	Yes
Amazon S3	-	-	Yes	Yes
Backblaze B2	-	-	Yes	Yes
Dropbox	-	-	Yes	Yes
Microsoft OneDrive	-	-	Yes	Yes
Google Drive	-	-	Yes	Yes
Google Cloud Storage	-	-	Yes	Yes
Microsoft Azure Blob + File Storage	-	-	Yes	Yes
WebDAV	-	-	Yes	Yes
OpenStack Swift	-	-	Yes	Yes
Box	-	-	Yes	Yes
Site Manager synchronization	-	-	Yes	Yes
Command-line interface	-	-	-	Yes
Batch transfers	-	-	-	Yes

Download Select Select Select

